# Chapter 6

## H200 Wireless Control Unit

#### The Control Unit is used to:

- Turn on/off the NESS H200 Wireless System.
- Test stimulation in the Orthosis.
- Select a stimulation program.
- Turn on/off and pause stimulation.
- Adjust the stimulation intensity level.
- Monitor system status.
- Mute system audio alerts.

The Control Unit is powered by a single rechargeable NiMH AAA battery.



CAUTION: Use only a battery supplied by Bioness Inc.

## H200° Wireless

### **Operating Buttons**

The Control Unit has eight operating buttons. See Table 6-1, Figure 6-1, and Table 6-2.

Operating Button	Description	Function
	On/Off Button	Turns On/Off the System
(	Trigger Button	Turns On/Off and Pauses Stimulation
+	Intensity Adjustment Buttons	Increases Stimulation Intensity Decreases Stimulation Intensity
	Mute Button	Mutes/Un-Mutes the Control Unit and H200 Wireless Orthosis Audio Alerts
[1 → <b>■</b> [2 → <b>■</b>	Program Selection Buttons	Top: Selects Program 1 Bottom: Selects Program 2
	Stimulation Test Button	Tests Stimulation in the H200 Wireless Orthosis

Table 6-1: Operating buttons: functions.





Figure 6-1: Operating buttons.

Visual Display	Description	Definition
	On/Off Button FLASHES GREEN	System On
	Trigger Button FLASHES YELLOW RAPIDLY	Stimulation On
(*************************************	Trigger Button is SOLID YELLOW	Stimulation Paused

Table 6-2: Operating buttons: visual displays.

## H200° Wireless

## **Indicators and Digital Display**

#### The Control Unit front panel features:

- A Control Unit status indicator: (i)
- An RF communication indicator:
- Program 1 and Program 2 selection indicators: 1
- A digital display. See Figure 6-2.



Figure 6-2: H200 Wireless Control Unit digital display and indicators.

![](_page_4_Picture_0.jpeg)

#### The Control Unit indicators show:

- Control Unit status.
- Stimulation program selected.
- Control Unit low battery.
- RF communication status. See Table 6-3.

Indicator	Description	Definition
	Control Unit Status Indicator FLASHES YELLOW	Low Battery, H200 Wireless Control Unit
i	Control Unit Status Indicator is SOLID RED	Control Unit Charging Error; Electronic Registration Error; Control Unit Hardware/Software Error
[1	Program 1 Indicator is GREEN	Program 1 Selected
2	Program 2 Indicator is GREEN	Program 2 Selected
	RF Communication Indicator FLASHES RED	RF Communication Error

Table 6-3: H200 Wireless Control Unit indicator lights.

## H200° Wireless

#### The Control Unit digital display shows:

- Stimulation intensity level: 0–9.
- Electronic registration status.
- Battery charging status. See Table 6-4.

Display	Description	Definition
	0–9	Stimulation Intensity Level; "0" Equals No Stimulation
Ī	Alternating GREEN Arches	Registration in Process
	Letter "C"	Registration Complete
E	Letter "E"	Registration Error
	Letter "U"	H200 Wireless Control Unit Unregistered
Ĉ	Rotating GREEN Circle	H200 Wireless Control Unit Charging
-	Horizontal GREEN Line	H200 Wireless Control Unit Fully Charged

Table 6-4: H200 Wireless Control Unit digital displays.

![](_page_6_Picture_0.jpeg)

#### **Audio Alerts**

#### The H200 Wireless Control Unit beeps to indicate:

- Electronic registration was initiated, successful, or unsuccessful.
- The NESS H200 Wireless System was turned on/off.
- A program has ended and stimulation has stopped.
- An Control Unit hardware/software error occurred.
- The Control Unit battery charge level is low.
- The audio alerts were muted/un-muted.
- A charger was connected.
- RF communication failed.
- A button was pressed.

## **Charging Port**

The Control Unit charging port is located at the bottom of the Control Unit, under the flexible cover. See Figure 6-3.

**Note:** Next to the charging port is the signal input/output port for the clinician's programmer.

![](_page_7_Picture_4.jpeg)

Figure 6-3: Control Unit charging port.

![](_page_7_Picture_6.jpeg)

CAUTION: The signal input/output port on the Control Unit is only to be used by the clinician during setup.