# Thigh Plus/Stand-Alone User's Reference Card



#### **Putting On the L300 Go Thigh Cuff**

- Wet and attach the Thigh Cloth Electrodes to the L300 Go Thigh Cuff.
- 2 For placement on the quadriceps while seated, place the L300 Go Thigh Cuff on your thigh, positioning the locator three finger widths from your knee.
- 3 Make sure to place the L300 Go Thigh Cuff in the fitting position determined by your clinician. See Figure 3.
- For placement on hamstrings, while seated or standing, place the L300 Go Thigh Cuff at least four finger widths above back of knee, right under the gluteal fold. See Figure 4.










### **Operating the L300 Go System**



- 1 Press the on/off button to turn on the EPG.
- Select either gait or training mode:
  - To select gait mode, press the mode button briefly.
  - To select training mode, press and *hold* the mode button for 3 seconds.

    The EPG will acknowledge the new mode with audio and tactile feedback. To exit training mode, press and hold the mode button for 3 seconds.
- 3 If necessary, adjust the stimulation intensity level by pressing the + plus or minus button.


#### Rx Only (US only)

For a complete list of all contraindications, warnings, adverse reactions, precautions, and care instructions please review the User's Guide.



#### **Charging the Batteries**

Connect the system charger set to the EPG, then plug into a wall socket.



#### **EPG Displays**



# **Control Unit Displays\***



	Solid	System is On
2	Stim Wave Visible	System is in Gait/Training, Stimulation is enabled but Off
2	Stim Wave Flashes	System is in Gait/Training, Stimulation is enabled and On
048	Displays 0-9	Intensity Level
ķ	Displays Gait Icon	Gait Mode is Enabled
<b>(H)</b>	Displays Training Icon	Training Mode is Enabled
<b>-</b>	Displays Audio Icon	Audio is On

<sup>\*</sup>Optional Component

## **System Error Displays**

